

Recombinant Mouse C-C Motif Chemokine 9,CCL9,,MIP-1-γ

Catalog No: #AP75912

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75912-1 10ug #AP75912-2 50ug #AP75912-3 500ug #AP75912-4 1mg

Description

Product Name	Recombinant Mouse C-C Motif Chemokine 9,CCL9,,MIP-1-γ
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P51670
Uniprot	P51670
GeneID	20308;
Calculated MW	11.6kDa
Target Sequence	MQITHATETKEVQSSLKAQQGLEIEMFHMGMFQDSSDCCLSYNSRIQCSRFYFPTSGGCTRPGIIFISKRGFQ VCANPSDRRVQRCIERLEQNSQPRTYKQ
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM TrisHCl300mM NaCl,pH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only