

# Recombinant Human UFM1-Activating Enzyme 5,UBA5 (N-6His)

Catalog No: #AP75913

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75913-1 10ug #AP75913-2 50ug #AP75913-3 500ug #AP75913-4 1mg

## Description

Product Name	Recombinant Human UFM1-Activating Enzyme 5,UBA5 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9GZZ9
Uniprot	Q9GZZ9
GeneID	79876;
Calculated MW	47kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMAESVERLQQRVQELERELAQERSLQVPRSGDGGGGRVRIEKMSSEVVDS NPYSRLMALKRMGIVSDYEKIRTFVAIVGVGGVGSVTAEMLTRCGIGKLLLFDYDKVELANMNRLLFFQPHQA GLSKVQAAEHTLRNINPDVLFVHNYNITTVENFQHFMDRISNGGLEEGKPVLDVLSVDNFEARMTINTACNE LGQTWMEGVSSENAVSGHIQLIIPGESACFACAPPLVVAANIDEKTLKREGVCAASLPTTMGVVAGILVQNVLK FLLNFGTVSFYLYGNAMQDFPTMSMKPNPQCDDRNCRKQQEYKKAALPKQEVQEEEEIIHEDNEWGIE LVSEVSEELKNFSGPVPDLPEGITVAYTIPKKQEDSVTELTVEDSGESLEDLMAKMKNM
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl pH 8.0, 1mM DTT, 10% glycerol, 50mM NaCl.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only