

Recombinant Human Small Ubiquitin-Related Modifier 3,SUMO3,SMT3A



Catalog No: #AP75916

Orders: order@signalwayantibody.com

Package Size: #AP75916-1 10ug #AP75916-2 50ug #AP75916-3 500ug #AP75916-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Small Ubiquitin-Related Modifier 3,SUMO3,SMT3A
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P55854
Uniprot	P55854
GeneID	6612;
Calculated MW	11.6kDa
Target Sequence	MSEEKPKEGVKTENDHINLKVAGQDGSVVQFKIKRHTPLSKLMKAYCERQGLSMRQIRFRFDGQPINETDTPA QLEMEDEDTIDVFQQQTGGVPESLAGHSF
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only