

Recombinant Soybean Birkproteinase Inhibitor,BBI (C-6His)

Catalog No: #AP75919

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75919-1 10ug #AP75919-2 50ug #AP75919-3 500ug #AP75919-4 1mg

Description

Product Name	Recombinant Soybean Birkproteinase Inhibitor,BBI (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Soybean
Accession No.	P01055
Uniprot	P01055
GeneID	548083;
Calculated MW	9kDa
Target Sequence	MDDESSKPCCDQCACTKSNPPQCRCSDMRLNSCHSACKSCICALSYPAQCFCVDITDFCYEPCKPSEDDKE NLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only