

Recombinant Human Disks Large Homolog 4,DLG4,PSD95 (N-6His)



Catalog No: #AP75922

Orders: order@signalwayantibody.com

Package Size: #AP75922-1 10ug #AP75922-2 50ug #AP75922-3 500ug #AP75922-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Disks Large Homolog 4,DLG4,PSD95 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P78352
Uniprot	P78352
GeneID	1742;
Calculated MW	81.9kDa
Target Sequence	MNHKVVHHHHHMDCLCIVTTKKYRYQDEDTPPLEHSPAHLPNQANSPPVIVNTDTLEAPGYELQVNGTEGEM EYEEITLERNGLGFSIAGGTDNPHIGDDPSIFITKIIPGGAAAQDGRRLRVNDSILFVNEVDVREVTHSAAVEAL KEAGSIVRLYVMRRKPPAEKVMEIKLIKGPKGLGFSIAGGVGNQHIPGDNSIYVTKIIEGGAAHKDGRLLQIGDKIL AVNSVGLLEDVMHEDAVAALKNTYDVVYLKVAKPSNAYLSDSYAPPDITTSYQHLNDNEISHSSYLGTDYPTAM TPTSPRRYSPVAKDLLGEEDIPREPRRIVHRGSTGLGFNIVGGEDGEGIFISFILAGGPADLSGELRKGDKILSV NGVDLRNASHEQAAIALKNAGQTVTIIAQYKPEEYSRFEAKIHDREQLMNSSLGSGTASLRSNPKRGFYIRAL FDYDKTKDCGFLSQALSFRFGDVLHVIDAGDEEWWQARRVHSDSETDDIGFIPSKRRVERREWSRLKAKDW GSSSGSQGREDSVLSYETVTQMEVHYARPIIILGPTKDRANDLLSEFPDKFGSCVPHTRPKREYEIDGRDY HFVSSREKMEKDIQAHKFIAGQYNSHLYGTSVQSVREVAEQGKHCILDVSANAVRRLQAAHLHPAIFIRPRS LENVLEINKRITEEQARKAFDRATKLEQEFTECFSAIVEGDSFEEIYHKVKRVIEDLSGPYIWVPARERL
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only