

Recombinant Mouse Adiponectin,Acrp30,AdipoQ ((N-6His, E. coli)

Catalog No: #AP75923

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75923-1 10ug #AP75923-2 50ug #AP75923-3 500ug #AP75923-4 1mg

Description

Product Name	Recombinant Mouse Adiponectin,Acrp30,AdipoQ ((N-6His, E. coli)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q60994
Uniprot	Q60994
GeneID	11450;
Calculated MW	27.2kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMEDDVTTEELAPALVPPKGTGAGWWMAGIPGHPGHNGTPGRDGRDGTGPEKGEKGDAGLLGPKGETGDVGMTGAEGPRGFPGTPGRKGEPGEAAVYRSAFSVGLETRVTVPNVPIRFTKIFYNQQNHYDYGSTGKFYCNIPGLYYFSYHITVYMKDVKVSFLFKKDKAVLFTYDQYQEKNVDAQSGSVLLHLEVG DQVWLQVYGDGDHNGLYADNVNDSTFTGFLYHDTN
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only