

Recombinant Human Tripartite Motif-Containing Protein 5, TRIM5, RNF88 (N-6His)

Catalog No: #AP75934

Orders: order@signalwayantibody.com

Package Size: #AP75934-1 10ug #AP75934-2 50ug #AP75934-3 500ug #AP75934-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Tripartite Motif-Containing Protein 5, TRIM5, RNF88 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9C035
Uniprot	Q9C035
GeneID	85363;
Calculated MW	30.8kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMASGILVNVKEEVTCPICLELLTQPLSLDCGHSFCQACL TANHKKSM LDKGE SSCPVCRISYQPENIRPNRHVANIVEKLREVKLSPEGQKVDHDCARHGKLLLFQEDGKVICWLCERSQEHRG HHTFLTEEVAREYQVKLQA ALEMLRQKQQAEELEADIREEKASWKTQIQYDKTNVLADFEQLRDILDWEESN ELQNLEKEEEDILKSLTNSETEMVQQTQSLRELISDLEHRLQGSMELLQ
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only