

Recombinant Human Omentin,Intelectin-1,ITLN-1 (N-6His)

Catalog No: #AP75936

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75936-1 10ug #AP75936-2 50ug #AP75936-3 500ug #AP75936-4 1mg

Description

Product Name	Recombinant Human Omentin,Intelectin-1,ITLN-1 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8WWA0
Uniprot	Q8WWA0
GeneID	55600;
Calculated MW	32.7kDa
Target Sequence	MNHKVVHHHHHMTDEANTYFKEWTCSSSPSLPRSCKEIKDECPSAFDGLYFLRTENGVIIYQTFCDMTSGGG GWTLVASVHENDMRGKCTVGDRWSSQQGSKAVYPEGDGNWANYNTFGSAEAATSDDYKNPGYYDIQAKDL GIWHVPNKSPMQHWRNSSLLRYRTDTGFLQTLGHNLFGIYQKYPVKYGEKGCWTDNGPVIPIVYDFGDAQK TASYSPYQGREFTAGFVQFRVFNNERAANALCAGMRVTGCNTEHHCIGGGGYFPEASPPQCCGDFSGFDW SGYGTHVGYS
Formulation	Lyophilized from a 0.2 µm filtered solution of 50mM Tris,100mMNaCl, 5mM GSH, 0.5mM GSSG,pH8.0 .
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only