

Recombinant Human Repulsive Guidance Molecule A, RGMA (N-6His)

Catalog No: #AP75937

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75937-1 10ug #AP75937-2 50ug #AP75937-3 500ug #AP75937-4 1mg

Description

Product Name	Recombinant Human Repulsive Guidance Molecule A, RGMA (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96B86
Uniprot	Q96B86
GeneID	56963;
Calculated MW	29.8kDa
Target Sequence	MNHKVVHHHHHMPHLRFTDRFQTCKVQGAWPLIDNNYLVNVTNTPVLPGSAATATSKLTIIFKNFQECVDQ KVYQAEMDELPAAFVDGSKNGGDKHGANSKITEKVSQGHVEIQAKYIGTTIVVRQVGRYLTFAVRMPEEVVN AVEDWDSQGLYLCLRGCPNQQIDFQAFHTNAEGTGARRLAAASPAPTAPETFPYETAVAKCKEKLPEVDLY YQACVFDLLTTGDVNFVLAAYYALDVKMLHSNKDKLHLYERTRDLPG
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only