

## Recombinant Human Neurtin, NRN1 (N-6His)

Catalog No: #AP75939



Package Size: #AP75939-1 10ug #AP75939-2 50ug #AP75939-3 500ug #AP75939-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Neurtin, NRN1 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9NPD7
Uniprot	Q9NPD7
GeneID	51299;
Calculated MW	12.1kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMAGKCDVFKGFSDCLLKLGDSMANYPQGLDDKTNIKTVCTYWEDFHSTV TALTDQCQEGAKDMWDKLRKESKNLNIQGSLELCGSGNG
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only