

Recombinant Human Ezrin,EZR,VIL2

Catalog No: #AP75940



Package Size: #AP75940-1 10ug #AP75940-2 50ug #AP75940-3 500ug #AP75940-4 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Ezrin,EZR,VIL2
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P15311
Uniprot	P15311
GeneID	7430;
Calculated MW	69.4kDa
Target Sequence	MPKPINVRVTMDAELEFAIQPNNTGKQLFDQVVKTIGLREVWYFGLHYVDNKGFPWLKLDKKVSAQEVKRNPLQFKFRAKFYPEDVAEELIQDITQKLFLLQVKEGILSDEIYCPPETAVLLGSYAVQAKFGDYNKEVHKSGYLSERLIPQRVMDQHKLTRDQWEDRIQVWHAEHRGMLKDNAMLEYLKIAQDLEMYGINYFEIKNKKGTDLWLGV DALGLNIYEKDDKLTTPKIGFPWSEIRNISFNDKFKVPIKPIDKKAPDFVFYAPRLRINKRILQLCMGNHEL YMRRRKPDTIEVQQMKAQAREEKHQKQLERQQLETEKKRRETVEREKEQMMREKEELMLRLQDYEEKTKKAERELSE QIQRALQLEEEERKRAQEEAERLEADRMAALRAKEELERQAVDQIKSQEQLAAELA EYTAKIALLEEARRKED EVEEWQHRAKEAQDDLVKTKHEELHLVMTAPPPPPPPVYEPVSYHVQESLQDEGAEPTGYS AELSSEGIRDDR NEEKRITEAEKNERVQRQLLTLSS ELSQARDENKRTHNDIIHNENMRQGRDKYKTLRQIRQGN TKQRIDEFEAL
Formulation	Lyophilized from a 0.2 µm filtered solution of 10mM Hepes, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only