

Recombinant Human Complement Component C8 γ Chain, C8G (N-6His)



Catalog No: #AP75941

Orders: order@signalwayantibody.com

Package Size: #AP75941-1 10ug #AP75941-2 50ug #AP75941-3 500ug #AP75941-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Complement Component C8 γ Chain, C8G (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P07360
Uniprot	P07360
GeneID	733;
Calculated MW	22.6kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMQKPQRPRRPASPISTIQPKANFDAQQFAGTWLLVAVGSACRFLQEQQHRA EATTLHVAPQGTAMAVSTFRKLDGICWQVRQLYGDTGVLGRFLLQARGARGAVNWWVAETDYQSFAVLYLER AGQLSVKLYARSLPVSDSVLSGFQQRVQEAHLTEDQIFYFPKYGFCEAADQFHVLDEVRR
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only