

# Recombinant Human Dipeptidyl-Peptidase 3,DPP3 (N-6His)



Catalog No: #AP75944

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75944-1 10ug #AP75944-2 50ug #AP75944-3 500ug #AP75944-4 1mg

## Description

Product Name	Recombinant Human Dipeptidyl-Peptidase 3,DPP3 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9NY33
Uniprot	Q9NY33
GeneID	10072;
Calculated MW	84.8kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMADTQYILPNDIGVSSLDCREAFRLSPETERLYAYHLSRAAWYGGGLAVLLQT SPEAPYIYALLSRLFRAQDPDQLRQHALLAEGLTEEEYQAFVLYAAGVYSNMGNYSFGDTKFVPLPKEKLER VILGSEAAQQHPEEVRLWQTCGELMFSLEPRLRHLGLGKEGITYFSGNCTMEDAKLAQDFLDSQNL SAYN TRLFKEVDGEGKPYEVR LASVLGSEPSLDSEVT SKLKS YEFRGSPFQVTRGDYAPILQKVVEQLEKAKAYAA NSHQGQMLAQYIESFTQGSIEAHKRGSRFWIQDKGPIVESYIGFIESYRDPFGSRGEFEGFVAVVNKAMSAKF ERLVASAEQLLKELPWPPTFEKDKFLTPDFTSLDVLTFAGSGIPAGINIPNYDDL RQTEGFKNVSLGNVLAVAYA TQREKLTFL EEDDKDLYILWKGPSFDVQVGLHELLGHGSGKLFVQDEKGFANFDQETVINPETGEQIQSWYRS GETWDSKFS TIASSYEECRAESVGLYLCLHPQVLEIFGFEGADAEDVIYVNWLN MVRAGLLALEFYTP EAFNW RQAHMQARFVILRVLLEAGEGLVTITPTTGSDGRPDARVRLDRSKIRSVGK PALERFLRRLQVLKSTGDVAGG RALYEGYATVTDAPPECFLTRD TVLLRKESRKLIVQPNTRLEGS DVQLLEYEASAAGLIRSF SERFPEDGP EL EEILTQLATADARFWKGPSEAPSGQA
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,pH8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only