

## Recombinant Human Ferritin Light Chain,FTL (N-6His)

Catalog No: #AP75946



Package Size: #AP75946-1 10ug #AP75946-2 50ug #AP75946-3 500ug #AP75946-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Ferritin Light Chain,FTL (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P02792
Uniprot	P02792
GeneID	2512;
Calculated MW	24.5kDa
Target Sequence	MNHKVVHHHHHMSQIRQNYSTDVEAAVNSLVNLYLQASYTYLSLGFYFDRDDVALEGVSHFFRELAEEKRE GYERLLKMQNQRRGGRALFQDIKPAEDEWGTKTPDAMKAAMALEKKLNQALLDLHALGSARTDPHLCDFLETH FLDEEVKLIKMGDHLTNLHRLGGPEAGLGEYLFERLTLKHD
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only