

Recombinant Human Interleukin-22,IL-22

Catalog No: #AP75947



Package Size: #AP75947-1 10ug #AP75947-2 50ug #AP75947-3 500ug #AP75947-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Interleukin-22,IL-22
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9GZX6
Uniprot	Q9GZX6
GeneID	50616;
Calculated MW	16.9kDa
Target Sequence	MAPISSHCRLDKSNFQQPYITNRTFMLAKEASLADNNTDVRLLIGEKLFGVSMSERCYLMKQVLNFTLEEVLP QSDRFQPYMQEVVPFLARLSNRLSTCHIEGDDLHIQRNVQKLKDTVKKLGESGEIKAIGELDLLFMSLRNACI
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only