

Recombinant Human Ferritin Heavy Chain,FTH (N-6His)

Catalog No: #AP75948

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75948-1 10ug #AP75948-2 50ug #AP75948-3 500ug #AP75948-4 1mg

Description

Product Name	Recombinant Human Ferritin Heavy Chain,FTH (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P02794
Uniprot	P02794
GeneID	2495;
Calculated MW	23.4kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMTTASTSQVRQNYHQDSEAAINRQINLELYASYVYLSMSYYFDRDDVALKNF AKYFLHQSHEEREHAEKLMKLNQRGGRIFLQDIKKPDCDDWESGLNAMECALHLEKNVNQSLLELHKLATD KNDPHLCDFIETHYLNEQVKAIKELGDHVTNLRKMGAPESGLAEYLFDKHTLGSDSNES
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only