

Recombinant Human Polyubiquitin-B,UBB (N-GST, 6His)

Catalog No: #AP75950

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75950-1 10ug #AP75950-2 50ug #AP75950-3 500ug #AP75950-4 1mg

Description

Product Name	Recombinant Human Polyubiquitin-B,UBB (N-GST, 6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P0CG47
Uniprot	P0CG47
GeneID	7314;
Calculated MW	40.9kDa
Target Sequence	MSPILGYWKIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWRNKKFELGLEFPNLPYYIDGDVKLTQSMARIYI ADKHNMLGGCPKERAIEISMLEGAVLDIRYGVSRAYSKDFETLKVDFLSKLPPEMLKMFEDRLCHKTYLNGDHV THPDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAIQIDKYLKSSKYIAWPLQGWQATFGGGDHPKSDG STSGSGHHHHHSAGLVPRGSTAIGMKETAALKFERQHMDSPDLGTGGGSGDDDDKSPMGYRGSMQIFVK TLTGKTITLEVEPSDTIENVKAKIQDKEGIPPDQQLIFAGKQLEDGRTLSDYNIQKESTLHLVLRRLGG
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only