

Recombinant Human Ubiquitin-Conjugating Enzyme E2 D1,UBE2D1,UbcH5a (N-GST)



Catalog No: #AP75951

Orders: order@signalwayantibody.com

Package Size: #AP75951-1 10ug #AP75951-2 50ug #AP75951-3 500ug #AP75951-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Ubiquitin-Conjugating Enzyme E2 D1,UBE2D1,UbcH5a (N-GST)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P51668
Uniprot	P51668
GeneID	7321;
Calculated MW	42.9kDa
Target Sequence	MSPILGYWKIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWRNKKFELGLEFPNLPYYIDGDVKLTQSMAIIRYI ADKHNMMLGGCPKERAEISMLEGAVLDIRYGVSRIAYSDFETLKVDFLSKLPPEMLKMFEDRLCHKTYLNGDHV THPDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAIQIDKYLKSSKYIAWPLQGWQATFGGGDHPKSDL VPRGSMALKRIQKELSDLQRDPPAHCSAGPVGDDLFHWQATIMGPPDSAYQGGVFFLTVHFPTDYPFKPPKI AFTTKIYHPNINSNGSICLDILRSQWSPALTVSKVLLSICSLLCDPNPDDPLVPDIAQIYKSDKEKYNRHAREWTQ KYAM
Formulation	Lyophilized from a 0.2 µm filtered solution of 50mM HEPES150mM NaCl2mM DTT,10% Glycerin pH7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only