

Recombinant Human Annexin A3,ANXA3

Catalog No: #AP75952



Package Size: #AP75952-1 10ug #AP75952-2 50ug #AP75952-3 500ug #AP75952-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Annexin A3,ANXA3
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P12429
Uniprot	P12429
GeneID	306;
Calculated MW	36.4kDa
Target Sequence	<p>MASIWVGHRTVDRDYPDFSPSVDAAEAIQKAIKIGTDEKMLISILTERSNAQRQLVKEYQAAYGKELKDDLKGDLSGHFEHLMVALVTPPAVFDKQLKSMKGAGTNEDALIEILTTRTSRQMKDISQAYTYVYKSLGDDISSETSGDFRKALLTLADGRRDESLKVDEHLAKQDAQILYKAGENRWGTDEDKFEILCLRSFPQLKLTDFEYRNISQKDIVDSIKGELSGHFEDLLLAINCVRNTPAFLAERLHRALKGIGTDEFTLNRIMVSRSEIDLLDIRTEFKKHGYSLYSAIKSDTSGDYEITLLKICGGDD</p>
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only