

Recombinant Human Peroxiredoxin-1, PRDX1 (N, C-6His)

Catalog No: #AP75953

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75953-1 10ug #AP75953-2 50ug #AP75953-3 500ug #AP75953-4 1mg

Description

Product Name	Recombinant Human Peroxiredoxin-1, PRDX1 (N, C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q06830
Uniprot	Q06830
GeneID	5052;
Calculated MW	25.3kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMSSGNAKIGHPAPNFKATAVMPDGQFKDISLSDYKGKYVFFFYPLDFTFVC PTEIIAFSDRAEEFKKLNQCQVIGASVDSHFCHLAWVNTPKKQGGLGPMNIPLVSDPKRTIAQDYGVKKADEGIS FRGLFIIDDKGILRQITVNDLPVGRSVDLRLVQAFQFTDKHGEVCPAGWKPGSDTIKPDVQKSKEYFSKQKL EHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of PBS, 10% glycerol, 0.1mM DTT, pH 6.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only