

Recombinant Human Transcriptional Repressor Protein YY1 (C-6His)

Catalog No: #AP75954

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75954-1 10ug #AP75954-2 50ug #AP75954-3 500ug #AP75954-4 1mg

Description

Product Name	Recombinant Human Transcriptional Repressor Protein YY1 (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P25490
Uniprot	P25490
GeneID	7528;
Calculated MW	12.6kDa
Target Sequence	MVTMWSSDEKKDIDHETVVEEQIIGENSPDYSEYMTGKKLPPGGIPGIDLSDPKQLAEFARMKPRKIKEDDAP RTIACPHKGCTKMFRDNSAMRKHLHHTHGLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only