

Recombinant Human Peptidyl-Prolyl Cis-Trans Isomerase FKBP3 (N-6His)



Catalog No: #AP75955

Orders: order@signalwayantibody.com

Package Size: #AP75955-1 10ug #AP75955-2 50ug #AP75955-3 500ug #AP75955-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Peptidyl-Prolyl Cis-Trans Isomerase FKBP3 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q00688
Uniprot	Q00688
GeneID	2287;
Calculated MW	27.3kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMAAAVPQRAWTVELRSEQLPKKDIKFLQEHGSDSFLAEHKLLGNIKNVAK TANKDHLVTAYNHLFETKRFKGTESISKVSEQVKNVKNLNEKPKETKSEETLDEGPPKYTKSVLKKGDKTNFP KKGDVVHCWYTGTLQDGTVFDTNIQTSAKKKKNAKPLSFKVGVGKIVIRGWDEALLTMSKGEKARLEIEPEWA YGKKGQPDAKIPPNAKLTFEVELVDID
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl1mM DTT,10% glycerol .
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only