

# Recombinant Human Activating Transcription Factor 1, ATF1 (C-6His)

Catalog No: #AP75956

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75956-1 10ug #AP75956-2 50ug #AP75956-3 500ug #AP75956-4 1mg

## Description

Product Name	Recombinant Human Activating Transcription Factor 1, ATF1 (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P18846
Uniprot	P18846
GeneID	466;
Calculated MW	30.3kDa
Target Sequence	MEDSHKSTTSETAPQPGSAVQGAHISHIAQQVSSLSESEESQDSSDSIGSSQKAHGILARRPSYRKILKDLSSSE DTRGRKGDGENSGVSAAVTSMSVPTPIYQTSSGQYIAIAPNGALQLASPGTDGVQGLQTLTMTNSGSTQQGT TILQYAQTSQGQILVPSNQVVVQTASGDMQTYQIRTTPSATSLPQTVVMTSPVTLTSQTTKTDDPQLKREIRL MKNREAARECRRKKKEYVKLENRVAVLENQNKTLIEELKTLKDLYSNKSVVEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only