

# Recombinant Human Heat Shock Protein β-2,,HSPB2,MKBP (C-6His)

Catalog No: #AP75958

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75958-1 10ug #AP75958-2 50ug #AP75958-3 500ug #AP75958-4 1mg

## Description

Product Name	Recombinant Human Heat Shock Protein β-2,,HSPB2,MKBP (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q16082
Uniprot	Q16082
GeneID	3316;
Calculated MW	21.3kDa
Target Sequence	MSGRSVPHAHPATAEYEFANPSRLGEQRFGEGLLP EEILTPTLYHGYVVRPRAAPAGEGSRAGASELRLSEG KFQAFLDVSHFTPDEVTVRTVDNLLLEVSARHPQRLDRHG FVSREFCRTYVLPADVDPWRVRAALSHDGILNLE APRGGRHLDTEVNEVYISLLPAPPDPEEEEEAAIVEPLEHHHHHH
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM HEPES,0.1MKCl,pH7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only