

Recombinant Mouse Sonic Hedgehog,SHH (C25II)

Catalog No: #AP75960



Package Size: #AP75960-1 10ug #AP75960-2 50ug #AP75960-3 500ug #AP75960-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Sonic Hedgehog,SHH (C25II)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q62226
Uniprot	Q62226
GeneID	20423;
Calculated MW	19.8kDa
Target Sequence	MIIGPGRGFGKRRHPKLTPLAYKQFIPNVAEKT LGASGRYEGKITRNSERFKELTPNYNPDIIIFKDEENTGADR LMTQRCKDKLNALAISVMNQWPGVKLRVTEGWDEDGHHSEESLHYEGRAVDITTSRDRSKYGMLARLAVE AGFDWVYYESKAHIHCSVKAENSVAAKSGG
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 1mM DTT, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only