

## Recombinant Mouse Sonic Hedgehog, SHH

Catalog No: #AP75963



Package Size: #AP75963-1 10ug #AP75963-2 50ug #AP75963-3 500ug #AP75963-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse Sonic Hedgehog, SHH
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q62226
Uniprot	Q62226
GeneID	20423;
Calculated MW	19.8kDa
Target Sequence	MCGPGRGFGKRRHPKLLTPLAYKQFIPNVAEKT LGASGRYEGKITRNSERFKELTPNYNPDIIIFKDEENTGAD RLMTQRCKDKLNLALAISVMNQWPGVKLRVTEGWDEDGHHSEESLHYEGRAVDITTSRDRDRSKYGMLARLAV EAGFDWVYYESKAHIHCSVKAENSVAAKSGG
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 1mM DTT, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only