

## Recombinant Human Protachykinin-1, TAC1 (N-6His)

Catalog No: #AP75965



Package Size: #AP75965-1 10ug #AP75965-2 50ug #AP75965-3 500ug #AP75965-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Protachykinin-1, TAC1 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P20366
Uniprot	P20366
GeneID	6863;
Calculated MW	15.2kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMEEIGANDDLNYSDWYDSDQIKEELPEPFHELLQRIARRPKPQQFFGLMG KRDADSSIEKQVALLKALYGHGQISHKRHKTDSEVGLMGKRALNSVAYERSAMQNYERRR
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH 8.0 (10% glycerol 2M Urea).
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only