

# Recombinant Human Heat Shock Protein $\beta$ -8,,HSPB8 (C-6His)

Catalog No: #AP75967

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75967-1 10ug #AP75967-2 50ug #AP75967-3 500ug #AP75967-4 1mg

## Description

Product Name	Recombinant Human Heat Shock Protein $\beta$ -8,,HSPB8 (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UJY1
Uniprot	Q9UJY1
GeneID	26353;
Calculated MW	22.7kDa
Target Sequence	MADGQMPPFSCHYPSRLRRDPFRDSPLSSRLDDGFGMDPFPDDL TASWPDWALPRLSSAWPGTLRSGMVP RGPTATARFGVPAEGRTPPPPFGEPWKVCVNVHSFKPEELMVKTKDGYVEVSGKHEEKQQEGGIVSKNFTK KIQLPAEVDPVTVFASLSPEGLLIIEAPQVPPYSTFGESSFNNELPQDSQEVTCTLEHHHHHHH
Formulation	Supplied as a 0.2 $\mu$ m filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only