

# Recombinant Human Glia Maturation Factor β,GMF-β,GMFB (C-6His)

Catalog No: #AP75968

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75968-1 10ug #AP75968-2 50ug #AP75968-3 500ug #AP75968-4 1mg

## Description

Product Name	Recombinant Human Glia Maturation Factor β,GMF-β,GMFB (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P60983
Uniprot	P60983
GeneID	2764;
Calculated MW	17.7kDa
Target Sequence	MSESLVVCDAEDLVEKLRKFRFRKETNNAIIMKIDKDKRLVVLDEELEGISPDELKDELPERQPRFIVYSYKY QHDDGRVSYPLCFIFSSPVGCKPEQQMMYAGSKNKLVQTAELTKVFEIRNTEDLTEEWLREKLGFFHLEHHH HHH
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM Tris,200mMNaCl,pH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only