

# Recombinant Human IA2,PTPRN (Arg576-Gln950, N-6His)

Catalog No: #AP75971

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75971-1 10ug #AP75971-2 50ug #AP75971-3 500ug #AP75971-4 1mg

## Description

Product Name	Recombinant Human IA2,PTPRN (Arg576-Gln950, N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q16849
Uniprot	Q16849
GeneID	5798;
Calculated MW	44.6kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMRQQDKERLAALGPEGAHGDTTFEYQDLCRQHMTKSLFNRAEGPPEPSR VSSVSSQFSDAAQASPSHSSTPSWCEEPAQANMDISTGHMILAYMEDHLRNRDLAKEWQALCAYQAEPN TCATAQGEQNIKKNRHPDFLPYDHARIKLVESSPSRSDYINASPIIEHDPRMPAYIATQGPLSHTIADFWQMV WESGCTVIVMLTPLVEDGVKQCDRYWPDEGASLYHYEVNLVSEHIWCEDFLVRSFYLNKVNQTQETRTLQF HFLSWPAEGTPASTRPLLDFFRRKVNKCYRGRSCPIIVHCSDGAGRTGTYLIDMVLNRMAGVKEIDIAATLEH VRDQRPLVRSKDQFEFALTAVAAEEVNAILKALPQ
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,pH8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only