

# Recombinant Human Pyruvate Kinase M2,PKM2 (N-6His)

Catalog No: #AP75972

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75972-1 10ug #AP75972-2 50ug #AP75972-3 500ug #AP75972-4 1mg

## Description

Product Name	Recombinant Human Pyruvate Kinase M2,PKM2 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P14618
Uniprot	P14618
GeneID	5315;
Calculated MW	60.1kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMSKPHSEAGTAFIQTQQLHAAMADTFLEHMCRLDIDSPITARNTGIICTIGPA SRSVETLKEMIKSGMNVARLNFSHGTHEYHAETIKNVRTATESFASDPILYRPVAVALDTKGPEIRTLGIKGGST AEVELKKGATLKITLDNAYMEKCDENILWLDYKNICKVVEVGSKIYVDDGLISLQVKQKQKADFLVTEVENGSSL GSKKGVNLPGAAVDLPVSEKDIQDLKFGVEQDVMVFASFIRKASDVHEVRKVLGEKGNIKIISKIENHEGV RRFDEILEASDGIMVARGDLGIEIPAEKVFLAQMMIGRCNRAGKPVICATQMLESMIKKPRPTRAEGSDVANA VLDGADCIMLSGETAKGDYPLEAVRMQHIAIEAEAAIYHLQLFEELRRLAPITSDPTEATAVGAVEASFKCCS GAIIVLTKSGRSAHQVARYRPRAPIIIVTRNPQTARQAHLRYGIFPVLCQDPVQEAWAEDVLDLRVNFAMNVGK ARGFFKKGDVIVLTGWRPGSGFTNTMRVVPVP
Formulation	Supplied as a 0.2 µm filtered solution of PBS, pH7.0, 10% glycerol.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only