

Recombinant Human IA2,PTPRN (Gln607-Asn686 & Trp795-Leu888, N-6His)



Catalog No: #AP75973

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75973-1 10ug #AP75973-2 50ug #AP75973-3 500ug #AP75973-4 1mg

Description

Product Name	Recombinant Human IA2,PTPRN (Gln607-Asn686 & Trp795-Leu888, N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q16849
Uniprot	Q16849
GeneID	5798;
Calculated MW	22.9kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMQDKERLAALGPEGAHGDTTFEYQDLCRQHMATKSLFNRAEGPPEPSRVS SVSSQFSDAAQASPSSHSSTPSWCEEPAQANGGGGSGGGGSGGGGSWQMWWESGCTVIVMLTPLVEDGV KQCDRYWPDEGASLYHYEVNLVSEHIWCEDFLVRSFYLNKVNQTQETRRTLQFHFLSWPAEGTPASTRPLL
Formulation	Supplied as a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only