

Recombinant Human Pro-Neuregulin-1, NRG1 β 1, HRG1 β 1 (Thr176-Lys246)

Catalog No: #AP75974

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75974-1 10ug #AP75974-2 50ug #AP75974-3 500ug #AP75974-4 1mg

Description

Product Name	Recombinant Human Pro-Neuregulin-1, NRG1 β 1, HRG1 β 1 (Thr176-Lys246)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q02297-6
Uniprot	Q02297
GeneID	3084;
Calculated MW	8.2kDa
Target Sequence	TSHLVKCAEKEKTFVNGGECFMVKDLSNPSRYLCKCPNEFTGDRQCNYVMASFYKHLGIEFMEAEELYQK
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only