

Recombinant Human Fibroblast Growth Factor 2,FGF-2,bFGF (K128N)

Catalog No: #AP75976

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75976-1 10ug #AP75976-2 50ug #AP75976-3 500ug #AP75976-4 1mg

Description

Product Name	Recombinant Human Fibroblast Growth Factor 2,FGF-2,bFGF (K128N)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P09038-2
Uniprot	P09038
GeneID	2247;
Calculated MW	17.2kDa
Target Sequence	MAAGSITTLPALPEDGGSGAFPPGHFKDPKRLYCKNGGFFLRIHPDGRVDGVREKSDPHIKLQLQAEERGVS IKGVCANRYLAMKEDGRLLASKCVTDECFERLESNNYNTYRSRKYTSWYVALNRTGQYKLGSKTGPGQKA ILFLPMSAKS
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris,200mM NaCl, pH7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 10 months.

Note: This product is for in vitro research use only