

Recombinant Human Persephin

Catalog No: #AP75978



Package Size: #AP75978-1 10ug #AP75978-2 50ug #AP75978-3 500ug #AP75978-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Persephin
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O60542
Uniprot	O60542
GeneID	5623;
Calculated MW	10.4kDa
Target Sequence	GSALSGPCQLWSLTLVAELGLGYASEEKVIFRYCAGSCPRGARTQHGLALARLQGQGRAHGGPCCRPTRY TDVAFLDDRHRWQRLPQLSAAACGCGG
Formulation	Lyophilized from a 0.2 µm filtered solution of 4mM HCl.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 12 months.

Note: This product is for in vitro research use only