

Recombinant Cavia porcellus CTLA-4 (C-6His)

Catalog No: #AP75979



Package Size: #AP75979-1 10ug #AP75979-2 50ug #AP75979-3 500ug #AP75979-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------|---|
| Product Name | Recombinant Cavia porcellus CTLA-4 (C-6His) |
| Host Species | P.Pichia |
| Purification | Greater than 95% as determined by reducing SDS-PAGE. |
| Specificity | Cavia porcellus |
| Accession No. | H0VUB1 |
| Uniprot | H0VUB1 |
| Calculated MW | 14.1kDa |
| Target Sequence | AMHVAQPAVVLAASSRGVASFCEYASSHNANEVRVTVLQQVASRTTEICAATYTVRELAFPEDSACAGTSS GTRVNLTIQGLRAADTGLYICKVELMYPPPYFVGTGNGTQIYVIDPEPCPDSDDHHHHHH |
| Formulation | Supplied as a 0.2 µm filtered solution of PBS, pH7.4. |
| Storage | Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles. |

Note: This product is for in vitro research use only