

## Recombinant Human Cyclin-H,CCNH (N-6His)

Catalog No: #AP75980



Package Size: #AP75980-1 10ug #AP75980-2 50ug #AP75980-3 500ug #AP75980-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Cyclin-H,CCNH (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P51946
Uniprot	P51946
GeneID	902;
Calculated MW	39.8kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMYHNSSQKRHWTFSSSEQLARLRADANRKFRCCKAVANGKVLNDP VFLEP HEEMTLCKYYEKRLLEFCSVFKPAMP RSVVGTACMYFKRFYLNNSVMEYHPRIIMLTCAFLACKVDEFNVSSP QFVGNLRESPLGQEKALEQILEYELLIIQLNFHLIVHNPYRPFEGFLIDLKTRYP ILENPEILRKTADDFLNRIAL TDAYLLYTPSQIALTAILSSASRAGITMESYLSESLMLKENRTCLSQLLDIMKSMRNLVKKYEPPRSEEVA VLKQ KLERCHSAELALNVITKKRKG YEDDDYVSKKSKHEEEEWTDDDL VESL
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris,100mM NaCl,2mM EDTA,2mM DTT,30%Glycerol,pH8.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only