

Recombinant Human Acylphosphate Phosphohydrolase 1, ACYP1 (N-6His)

Catalog No: #AP75981

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75981-1 10ug #AP75981-2 50ug #AP75981-3 500ug #AP75981-4 1mg

Description

Product Name	Recombinant Human Acylphosphate Phosphohydrolase 1, ACYP1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P07311
Uniprot	P07311
GeneID	97;
Calculated MW	13.42kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMAEGNTLISVDYEIFGKVQGVFFRKHTQAEGKKLGLVGWVQNTDRGTVQG QLQGPISKVRHMQEWLETRGSPKSHIDKANFNNEKVLKLDYSDFQIVK
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris, 100mM NaCl, 20% Glycerol, 1mM DTT, pH8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only