

# Recombinant Human S100 Calcium Binding Protein A6,S100A6 (N-6His)



Catalog No: #AP75982

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75982-1 10ug #AP75982-2 50ug #AP75982-3 500ug #AP75982-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                 |   |
|-----------------|---|
| Product Name    | Recombinant Human S100 Calcium Binding Protein A6,S100A6 (N-6His)   |
| Host Species    | E.coli  |
| Purification    | Greater than 95% as determined by reducing SDS-PAGE.  |
| Specificity     | Human   |
| Accession No.   | P06703  |
| Uniprot         | P06703  |
| GeneID          | 6277;   |
| Calculated MW   | 12.5kDa   |
| Target Sequence | MGSSHHHHHSSGLVPRGSHMMACPLDQAIGLLVAIFHKYSGREGDKHTLSKKELKELIQKELTIGSKLQDAEI<br>ARLMEDLDRNKDQEVNFQEYVTLGALALIYNEALKG |
| Formulation     | Supplied as a 0.2 µm filtered solution of PBS,pH7.4.  |
| Storage         | Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.                            |

Note: This product is for in vitro research use only