

Recombinant Human N-Acyl-Aromatic-L-Amino Acid Amidohydrolase,ACY3 (N-6His)



Catalog No: #AP75983

Orders: order@signalwayantibody.com

Package Size: #AP75983-1 10ug #AP75983-2 50ug #AP75983-3 500ug #AP75983-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human N-Acyl-Aromatic-L-Amino Acid Amidohydrolase,ACY3 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96HD9
Uniprot	Q96HD9
GeneID	91703;
Calculated MW	37.4kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMCSLPVPREPLRRVAVTGGTHGNEMSGVYLARHWLHAPAE LQRASFSAVP VLANPAATSGCRRYVDHDLNRTFTSSFLNSRPTPDDPYEVTRARELNQLLGPKASGQAFDFVLDLHNTTANM GTCLIAKSSHEVFAMHLRHLQLQYPELSCQVFLYQRSGEESYNLDSVAKNGLGLELGPQPQGVLRADIFSR MRTLVA TVLDFIELFNQGTAFPAFEMEAYRPVGVVDFPRTEAGHLAGTVHPQLQDRDFQPLQPGAPIFQMFS GEDLLYEGESTVYPVFINEAAYYEKGVAFVQTEKFTFTVPAMPALTPAPSPAS
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris,100mM NaCl,1mM DTT,10%Glycerol,pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only