

Recombinant Human Mammalian Ependymin-Related Protein 1,EPDR1,UCC1 (C-6His)



Catalog No: #AP75986

Orders: order@signalwayantibody.com

Package Size: #AP75986-1 10ug #AP75986-2 50ug #AP75986-3 500ug #AP75986-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Mammalian Ependymin-Related Protein 1,EPDR1,UCC1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UM22
Uniprot	Q9UM22
GeneID	54749;
Calculated MW	22.6kDa
Target Sequence	APRPCAPQQWEGRQVMYQQSSGRNSRALLSYDGLNQRVRLDERKALIPCKRLFYIILLYKDGVMFQIDQA TKQCSKMTLTQPWDPLDIPQNSTFEDQYSIGGPQEQITVQEWSDRKSARSYETWIGIYTVKDCYPVQETFTIN YSVILSTRFFDIQLGIKDPSVFTPPSTCQMAQLEKMSVDCSDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only