

# Recombinant Human CD27 Ligand,TNFSF7,CD70 (C-Fc)

Catalog No: #AP75988

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75988-1 10ug #AP75988-2 50ug #AP75988-3 500ug #AP75988-4 1mg

## Description

Product Name	Recombinant Human CD27 Ligand,TNFSF7,CD70 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P32970
Uniprot	P32970
GeneID	970;
Calculated MW	43.6kDa
Target Sequence	<p>QQQLPLESLGWDVAELQLNHTGPQQDPRLYWQGGPALGRSFLHGPELDKGQLRIHRDGIYMMVHIQVTLAICS</p> <p>STTASRHHPTTLAVGICSPASRSISLLRLSFHQGCTIASQRLTPLARGDTLCTNLGTLLPSRNTDETFFGVQW</p> <p>VRPVDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFN</p> <p>WYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQV</p> <p>YTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTPPVLDSDGSFFLYSKLTVDKSRWQQG</p> <p>NVFSCSVMEALHNHYTQKSLSLSPGK</p>
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only