

## Recombinant Human Cathepsin A,CTSA (C-6His)

Catalog No: #AP75993



Package Size: #AP75993-1 10ug #AP75993-2 50ug #AP75993-3 500ug #AP75993-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Cathepsin A,CTSA (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P10619
Uniprot	P10619
GeneID	5476;
Calculated MW	52.2kDa
Target Sequence	APDQDEIQRLPGLAKQPSFRQYSGYLKSGSGSKHLHYWFVESQKDPENSPVVLWLNGGPGCSSLDGLLTEHG PFLVQPDGVTLEYNPYSWNLIANVLYLESPAGVGFSSYDDKFYATNDTEVAQSNFEALQDFRFLFPEYKNNKL FLTGESYAGIYIPTLAVLVMQDPSMNLQGLAVGNLSSYEQNDNSLVYFAYYHGLLGNRLWSSLQTHCCSQN KCNFYDNKDLECVTNLQEVARIVGNSGLNIYNLYAPCAGGVPSHFYKEDTVVVQDLGNIFTRLPLKRMWHQA LLRSGDKVRMDPPCTNTTAASYLNNPYVRKALNIPEQLPQWDMCNFLVNLQYRRLYRSMNSQYLKLLSSQK YQILLYNGDVIDMACNFMGDEWFDVSLNQNKMEVQRRPWLVKYGDSEQIAGFVKEFSHIAFLTIKAGAHMVPT DKPLAAFTMFSRFLNKQPYHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only