

Recombinant Mouse Testican-3,SPOCK3 (C-6His)

Catalog No: #AP75995



Package Size: #AP75995-1 10ug #AP75995-2 50ug #AP75995-3 500ug #AP75995-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Testican-3,SPOCK3 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q8BKV0-1
Uniprot	Q8BKV0
GeneID	72902;
Calculated MW	47.9kDa
Target Sequence	AAAVAVAGGRSDGGNFLDEKQWLTTISQYDKEVGQWNKFRDEVEDDYFRTWNP GKPFQALDPAKDPCLKT KCSRHKVCITQDAQTALCISHRRLTHSMKEVGGSHKQWRGLPSSTCKPCPIAYASPVCGSDGHSYSSQCKLE YQACVLGKQISIKCEGRCPSPDKSMNIGRNVKRACSDLEFREVANRLRDWFKALHESGSQNKKTKALLRPE RSRFDTSILPICKDSLGMFNRDLDTNYDLLLDQSELGSIYLDKNEQCTKAFFNSCDTYKDSLISNNEWCYCFQR QQDPPCHTELSNIQKRQGKLLGQYIPLCEDGYYKPTQCHGSGVCWCVD RYGNEVVGSRINGVADCAID FEISGDFASGDFREWTDDEGEEDDIMNDKDDIEDDEDEGDDDDDDGDVHDGYIVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only