

Recombinant Mouse Transglutaminase 2,TGM2 (C-6His)

Catalog No: #AP75997

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75997-1 10ug #AP75997-2 50ug #AP75997-3 500ug #AP75997-4 1mg

Description

| | |
|-----------------|---|
| Product Name | Recombinant Mouse Transglutaminase 2,TGM2 (C-6His) |
| Host Species | Human Cells |
| Purification | Greater than 95% as determined by reducing SDS-PAGE. |
| Specificity | Mouse |
| Accession No. | P21981 |
| Uniprot | P21981 |
| GeneID | 21817; |
| Calculated MW | 78kDa |
| Target Sequence | AEELLERCdleIQANGRDHHTADLCQEKLVLRRGQRFRLTLyFEGRGYEASVDSLTFGAVTGPDPSEEAGTK ARFSLSDNVEEGSWSASVLDQQDNVLSLQLCTPANAPIGLYRLSLEASTGYQGSSFVLGHFILLYNWCPADD VYLDSEEEERREYVLTQQGFIYQGSVKFIKSVPNWFGQFEDGILDTCMLLDMPKFLKNRSRDCSRRSSPIYV GRVVSAMVNCNDDQGVLLGRWDNNYGDGISMWIGSVDILRRWKEHGCGQVQYQGCWVFAAVACTVLRRC LGIPTRVVTNYSAHQNSNLLIEYFRNEFGELESNKSEMIWNFHCWVESWMTRPDLQPGYEGWQAIDPTPQ EKSEGTYCCGPVSVRAIKEGDLSTKYDAPFVFAEVNADVVDWIRQEDGSVLKSINRSLVVGQKISTKSVGRDD REDITHYKYPEGSPEEREVFTKANHLNKLAEKEETGVAMRIRVGDMSMGNDFDVFAHIGNDTSETRECRLL LCARTVSYNGVLGPECGTEDINLTDPYSENSIPLRILYKYSGLTESNLIKVRGLLIEPAANSYLLAERDLYLE NPEIKIRVLGEPKQNRKLVAEVSLKNPLSDPLYDCIFTVEGAGLTKEQKSVEVSDPVPAGDLVKARVDLFPDIDIG LHKLVVNFQCDKLSVKGYRNVIIIGPAVDHHHHHH |
| Formulation | Supplied as a 0.2 µm filtered solution of 20mM TrisHCl,150mM NaCl,20%Glycerol,pH7.5. |
| Storage | Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles. |

Note: This product is for in vitro research use only