

Recombinant Mouse Collagen α -1(III) Chain, COL3A1 (C-6His)



Catalog No: #AP75998

Orders: order@signalwayantibody.com

Package Size: #AP75998-1 10ug #AP75998-2 50ug #AP75998-3 500ug #AP75998-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Collagen α -1(III) Chain, COL3A1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P08121
Uniprot	P08121
GeneID	12825;
Calculated MW	96.6kDa
Target Sequence	<p>QFDSYDVKSGVGGMGGYPGAPPPGPPGSSGHPGSPSPGYQGPPGEPGQAGPAGPPGPPGALGP AGPAGKDGESRPRGERGLPGPPGIKGPAGMPGFPGMKGHRGFDGRNGEKGETGAPGLKGENLPGD NGAPGPMGPRGAPGERGRPLGAAGARGNDGARGSDGQP GPPGPGTAGFPGSPGAKGEVGPAGSPG SNGSPGQRGEPGQGHAGAQQPPGPPGNNNGSPGGKGE MGPAGIPGAPGLIGARGPPGAGTNGIPGTRG PSGEPGKNGAKGEPGARGERGEAGSPGIPGPKGEDGKDGSPGEPGANGLPGAAGERGPSGFRGPAGPNI PGEKGGPPGERGGPGPAGPRGVAGEPGRDGTGGPGIRGMPGSPGGPGNDGKPGPPGSQGESRPGPPG PSGPRGQPGVMGFPGPKGNDGAPGKNGERGGPGGGLPGPAGKNGETGPQGGPPGTGPAGDKGDSGPP GPQGLQGIPGTGGPPGENGKPGEPGPKGEVGPAGPAGGKGDGAPGERGPPGTAGIPGARGGAGPPGPE GGKGPAGPPGPPGASGSPGLQGMPGERGGPGSPGPKGEKGEPEGGAGADGVPKDGPRGPAGPIGPPGP AGQPGDKGEGGSPGLPGIAGPRGGPGERGEHPPGPAGFPGAPQNGEPGAKGERGAPGEKGE GPPG PAGPTGSSGPAGPPGPGQVKGERGSPGGPGTAGFPGRGLPGPPGNNGNP GPPGSPGAPGKDGPPGPA GNSGSPGNPGIAGPKGDAGQPGEKGGPPGAQGGSPGPLGIAGLTGARGLAGPPGMPGPRGSPGPGQIKG ESGKPGASGHNGERPPGPQGLPGQP GTAGEPGRDGNPGSDGQPGRDGSPPGGKDRGENGSPGAPGAP GHPGPPGPVGP SGKSGDRGETGPAGPSGAPGAPARGAPGPQGRDGETGERG SNGIKGHRGFPGNP GPPGSPGAAGHQGAIGSPGAPGRGPVGP HPPGKDGTS GHPGPIGPPGPRGNRGERGSEGSFGHPGQP GPPGPPGAPGCCGGGAAAIAGVGGEKSGGFSPYYGVDHHHHHH</p>
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM HAc-NaAc, 150mM NaCl, pH4.5.
Storage	Store at -20°C , stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only