

Recombinant Human Magnesium-Dependent Phosphatase 1,MDP1 (C-6His)

Catalog No: #AP76000

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76000-1 10ug #AP76000-2 50ug #AP76000-3 500ug #AP76000-4 1mg

Description

Product Name	Recombinant Human Magnesium-Dependent Phosphatase 1,MDP1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q86V88
Uniprot	Q86V88
GeneID	145553;
Calculated MW	21.1kDa
Target Sequence	MARLPKLAVFDLDYTLWPFWVDTHVDPPFHKSSDGTVRDRRGQDVRLYPEVPEVLKRLQSLGVPGAAASRT SEIEGANQLLELFDLFRYFVHREIYPGSKITHFERLQQKTGIPFSQMIFDERRNIVDVSKLGVTCIHIQNGMNL QTLSQGLETFKAQTGPLRSSLEESPFEAVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 10%Glycerol, pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only