

## Recombinant Human PRADC1,PAP21 (C-6His)

Catalog No: #AP76001



Package Size: #AP76001-1 10ug #AP76001-2 50ug #AP76001-3 500ug #AP76001-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                 |  |
|-----------------|--|
| Product Name    | Recombinant Human PRADC1,PAP21 (C-6His)  |
| Host Species    | Human Cells  |
| Purification    | Greater than 95% as determined by reducing SDS-PAGE.   |
| Specificity     | Human  |
| Accession No.   | Q9BSG0   |
| Uniprot         | Q9BSG0   |
| GeneID          | 84279;   |
| Calculated MW   | 20kDa  |
| Target Sequence | HGFRIHDYLYFQVLSPGDIRYIFTATPAKDFGGIFHTRYEQIHLVPAEPPEACGELSNGFFIQDQIALVERGGCS<br>FLSKTRVVQEHEGGRAVIISDNAVNDNSFYVEMIQDSTQRTADIPALFLLGRDGYMIRRSLEQHGLPWAIISIPVN<br>VTSIPTFELLQPPWTFWVDHHHHHH   |
| Formulation     | Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl, pH7.4.  |
| Storage         | Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months. |

Note: This product is for in vitro research use only