

# Recombinant Human V-Set and Transmembrane Domain-Containing 2A,VSTM2A (C-6His)

Catalog No: #AP76003

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76003-1 10ug #AP76003-2 50ug #AP76003-3 500ug #AP76003-4 1mg

## Description

Product Name	Recombinant Human V-Set and Transmembrane Domain-Containing 2A,VSTM2A (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8TAG5-2
Uniprot	Q8TAG5
GeneID	222008;
Calculated MW	25.5kDa
Target Sequence	SQAKFTEFPRNVTATEGQNVEMSCAFQSGSASVYLEIQWWFLRGPEDLDPGAEGAGAQVKLLPDRDPDSDG TKISTVKVQGNDISHKQLQISKVRKKDEGLYECRVTDANYGELQEHLKQAYLKVNSHARRMQAFEASPMWL QDMKPRKNVSAIPSSIHGSAHQRTHTSTSPQVAKIPKQSPQSGMETHFEPFILPLTNAPQKGQSYRVDRFM NGDFVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only