

# Recombinant Human Leukotriene A-4 Hydrolase,LTA4H (C-6His)

Catalog No: #AP76004

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76004-1 10ug #AP76004-2 50ug #AP76004-3 500ug #AP76004-4 1mg

## Description

Product Name	Recombinant Human Leukotriene A-4 Hydrolase,LTA4H (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P09960
Uniprot	P09960
GeneID	4048;
Calculated MW	70.3kDa
Target Sequence	MPEIVDTCSLSPASVCRTKHLHLRCSVDFTRRTLGTAAALTVQSQEDNLRSLVLDTKDLTIEKVVINGQEVKY ALGERQSYKGGSPMEISLPIALSKNQEIVIEISFETSPKSSALQWLTPEQTSQKEHPYLFSSQCQAIHCRAILPCQD TPSVKLTYYAEVSVPKELVALMSAIRDGETPDPEDPSRKIYKFIQKVPICPYLIALVVGALSRQIGPRTLWSEK EQVEKSAYEFSETESMLKIAEDLGGPYVWQYDLLVLPSPFYGGMENPCLTFVTPTLLAGDKSLSNVIAHEIS HSWTGNLVTNKTWDHFWLNEGHTVYLERHICGRFLGKFRHFNALGGWQELQNSVKTFGETHPFTKLVVDL TDIDPDVAYSSVPYEKGFALLFYLEQLLGGPEIFLGFLLKAYVEKFSYKSITDDWKDFLYSYFKDKVDVLNQVD WNAWLYSPGLPPIKPNYDMLTNACIALSQRWITAKEDDLNSFNATDLKDLSSHQLNEFLAQTQRAPLPLGHI KRMQEVYNFNAINNSEIRFRWLRLCIQSKWEDAIPALPKMATEQGRMKFTRPLFKDLAAFDKSHDQAVRITYQE HKASMHPVTAMLVKGDKLVVDVHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM MES,2mM EDTA,20%Glycerol,pH5.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only